



o2si is an ISO 17025, Guide 34 and 9001:2008 accredited company

VOA Vials for Water Analysis



You can purchase preserved vials for methods 501, 502, 503, 524, 601, 602, 603, 624, 8010, 8020, 8030, or 8260.

Need something different?
Call customer service so we can meet your specific needs!

prep smart®

Water analysis just got smarter!

Unpreserved samples must be analyzed within 24 hours of collection. Don't waste time and money fussing around with concentrations and preservative options. No matter which purgeable organic compounds you are analyzing, you can count on us to provide you with an appropriately preserved vial. Our products will increase your holding time up to **14 days!**

Meets EPA Preservation Protocol

This product is manufactured using 40 mL amber, glass vials with polytetrafluoroethylene (PTFE)-faced septa.

2030 Savage Road
Charleston, SC 29407
Toll Free: 866.272.0932
www.o2si.com

